

RNF5 Polyclonal Antibody

Cat No: HR1AP10746

For research use only

Overview

Product Name	RNF5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	RNF5 G16 NG2 RMA1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ring finger protein 5(RNF5) Homo sapiens The protein encoded by this gene contains a RING finger, which is a motif known to be involved in protein-protein interactions. This protein is a membrane-bound ubiquitin ligase. It can regulate cell motility by targeting paxillin ubiquitination and altering the distribution and localization of paxillin in cytoplasm and cell focal adhesions. [provided by RefSeq, Jul 2008],