

# LGAT1 Polyclonal Antibody

Cat No: HR1AP10454

For research use only

## Overview

Product Name	LGAT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	LPGAT1 FAM34A KIAA0205
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lysophosphatidylglycerol acyltransferase 1(LPGAT1) Homo sapiens This gene encodes a member of the lysophospholipid acyltransferase family. The encoded protein catalyzes the reacylation of lysophosphatidylglycerol to phosphatidylglycerol, a membrane phospholipid that is an important precursor for the synthesis of cardiolipin. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],