



APLF Polyclonal Antibody

Cat No: HR1AP10021

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | APLF Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 60kDa |
| GeneID?Human? | APLF C2orf13 PALF XIP1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | Aprataxin and PNK-like factor (EC 4.2.99.18) (Apurinic-apyrimidinic endonuclease APLF) (PNK and APTX-like FHA domain-containing protein) (XRCC1-interacting protein 1) |
| Background | aprataxin and PNKP like factor(APLF) Homo sapiens C2ORF13 is a component of the cellular response to chromosomal DNA single- and double-strand breaks (Iles et al., 2007 [PubMed 17353262]).[supplied by OMIM, Mar 2008], |