

Ini1 Polyclonal Antibody

Cat No: HR1AP9641

For research use only

Overview

Product Name	Ini1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	SMARCB1 BAF47 INI1 SNF5L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SMARCB1 BAF47 INI1 SNF5L1
Background	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily b, member 1(SMARCB1) Homo sapiens The protein encoded by this gene is part of a complex that relieves repressive chromatin structures, allowing the transcriptional machinery to access its targets more effectively. The encoded nuclear protein may also bind to and enhance the DNA joining activity of HIV-1 integrase. This gene has been found to be a tumor suppressor, and mutations in it have been associated with malignant rhabdoid tumors. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Dec 2015],