

Hamartin Polyclonal Antibody

Cat No: HR1AP9589

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | Hamartin Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 130kDa |
| GeneID?Human? | TSC1 KIAA0243 TSC |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | tuberous sclerosis 1 |
| Background | tuberous sclerosis 1(TSC1) Homo sapiens This gene encodes a growth inhibitory protein thought to play a role in the stabilization of tuberin. Mutations in this gene have been associated with tuberous sclerosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2009], |