

Bovine Serum Albumin Polyclonal Antibody

Cat No: HR1AP2294

For research use only

Overview

Product Name	Bovine Serum Albumin Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Cow
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	66kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Bovine serum albumin (BSA) is the main protein synthesized in the liver that constitutes over 50% of total serum proteins. BSA is often used as a protein standard and can also be used in cell culture media, forensics, and foodstuffs.