

GATA-2/3 (Acetyl-Lys336/304) Polyclonal Antibody

Cat No: HR1AP9559

For research use only

Overview

Product Name	GATA-2/3 (Acetyl-Lys336/304) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human:K336/304,Mouse:K336/304,Rat:K336/305
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	GATA2/GATA3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GATA2
Background	GATA binding protein 2(GATA2) Homo sapiens This gene encodes a member of the GATA family of zinc-finger transcription factors that are named for the consensus nucleotide sequence they bind in the promoter regions of target genes. The encoded protein plays an essential role in regulating transcription of genes involved in the development and proliferation of hematopoietic and endocrine cell lineages. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Mar 2009],