

IL18 Polyclonal Antibody

Cat No: HR1AP11983

For research use only

Overview

Product Name	IL18 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	IL18 IGIF IL1F4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	interleukin 18(IL18) Homo sapiens The protein encoded by this gene is a proinflammatory cytokine that augments natural killer cell activity in spleen cells, and stimulates interferon gamma production in T-helper type I cells. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Aug 2011],