

mStrawberry Polyclonal Antibody

Cat No: HR1AP2277

For research use only

Overview

Product Name	mStrawberry Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Species independent
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mStrawberry is a bright, monomeric red fluorescent protein. mStrawberry excitation and emission maxima are 574 and 596 nm, respectively.