

## Tetranectin Polyclonal Antibody

Cat No: HR1AP9853

For research use only

### Overview

Product Name	Tetranectin Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CLEC3B TNA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Tetranectin (TN) (C-type lectin domain family 3 member B) (Plasminogen kringle 4-binding protein)
Background	function:Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.,mass spectrometry: PubMed:10614823,online information:Tetranectin,similarity:Contains 1 C-type lectin domain.,subunit:Homotrimer.,tissue specificity:Found in plasma.,