

## CUX1 Polyclonal Antibody

Cat No: HR1AP10085

For research use only

### Overview

Product Name	CUX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	165kDa
GeneID?Human?	CUX1 CUTL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cut like homeobox 1(CUX1) Homo sapiens The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, Feb 2011],