

GLI-1 Polyclonal Antibody

Cat No: HR1AP9927

For research use only

Overview

Product Name	GLI-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	120kDa
GeneID?Human?	GLI1 GLI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Zinc finger protein GLI1 (Glioma-associated oncogene) (Oncogene GLI)
Background	<p>GLI family zinc finger 1(GLI1) Homo sapiens This gene encodes a member of the Kruppel family of zinc finger proteins. The encoded transcription factor is activated by the sonic hedgehog signal transduction cascade and regulates stem cell proliferation. The activity and nuclear localization of this protein is negatively regulated by p53 in an inhibitory loop. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009],</p>