

## GST-Tag Polyclonal Antibody

Cat No: HR1AP2256

For research use only

### Overview

Product Name	GST-Tag Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Species independent
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>GST (Glutathione S-Transferase) is a widely used protein tag encoded by the <i>Schistosoma japonicum</i>. GST provides both an easily detectable tag and a simple purification process with little effect on the biological function of the protein of interest. Antibodies to GST are useful for checking protein expression both in plaques and on Western blots as well as for immunoaffinity purification of proteins expressed in cells.</p>