

BFP Polyclonal Antibody

Cat No: HR1AP2279

For research use only

Overview

Product Name	BFP Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Species independent
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Blue Fluorescent Protein (BFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.