

Bcl-6 Polyclonal Antibody

Cat No: HR1AP10045

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | Bcl-6 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 90kDa |
| GeneID?Human? | BCL6 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | <p>B-cell CLL/lymphoma 6(BCL6) Homo sapiens The protein encoded by this gene is a zinc finger transcription factor and contains an N-terminal POZ domain. This protein acts as a sequence-specific repressor of transcription, and has been shown to modulate the transcription of STAT-dependent IL-4 responses of B cells. This protein can interact with a variety of POZ-containing proteins that function as transcription corepressors. This gene is found to be frequently translocated and hypermutated in diffuse large-cell lymphoma (DLCL), and may be involved in the pathogenesis of DLCL. Alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, Aug 2015].</p> |