

# Myc-tag Polyclonal Antibody

Cat No: HR1AP2186

For research use only

## Overview

Product Name	Myc-tag Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Species independent
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	Myc Tag
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	c Myc; c Myc Epitope Tag; c Myc Tag; Myc Epitope Tag; Myc proto-oncogene protein; Transcription factor p64
Background	Myc tag is a polypeptide protein tag derived from the c-myc gene product that can be added to a protein. It can help to separate recombinant, overexpressed protein from wild type protein expressed by the host organism. Myc Tag also can be used to isolate protein complexes with multiple subunits.Myc Tag Antibody can be used in various immunoassays, such as ELISA, Western blotting, immunoprecipitation, immunofluorescence, and more.