

# PCNA Polyclonal Antibody

Cat No: HR1AP9950

For research use only

## Overview

Product Name	PCNA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	PCNA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Proliferating cell nuclear antigen (PCNA) (Cyclin)
Background	<p>proliferating cell nuclear antigen(PCNA) Homo sapiens The protein encoded by this gene is found in the nucleus and is a cofactor of DNA polymerase delta. The encoded protein acts as a homotrimer and helps increase the processivity of leading strand synthesis during DNA replication. In response to DNA damage, this protein is ubiquitinated and is involved in the RAD6-dependent DNA repair pathway. Two transcript variants encoding the same protein have been found for this gene. Pseudogenes of this gene have been described on chromosome 4 and on the X chromosome. [provided by RefSeq, Jul 2008],</p>