

FoxO3A Polyclonal Antibody

Cat No: HR1AP9934

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | FoxO3A Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 71kDa |
| GeneID?Human? | FOXO3 FKHRL1 FOXO3A |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | Forkhead box protein O3 (AF6q21 protein) (Forkhead in rhabdomyosarcoma-like 1) |
| Background | forkhead box O3(FOXO3) Homo sapiens This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq, Jul 2008], |