

HA-tag Polyclonal Antibody

Cat No: HR1AP2188

For research use only

Overview

Product Name	HA-tag Polyclonal Antibody
Source	Rabbit
Applications	WB,IP,ELISA
Species Reactivity	Species independent
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	HA Tag
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HA epitope tag; hemagglutinin
Background	The HA tag (Y-P-Y-D-V-P-D-Y-A) is derived from an epitope (amino acids 98-106) of the influenza hemagglutinin (HA) protein. HA tag is used as a general epitope tag in expression vectors for many recombinant proteins. The presence of HA tag facilitates the detection, isolation and purification of recombinant fusion proteins.