

## V1RL3 Polyclonal Antibody

Cat No: HR1AP6862

For research use only

### Overview

Product Name	V1RL3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	VN1R3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VN1R3; V1RL3; FKSG46; Vomeronasal type-1 receptor 3; V1r-like receptor 3
Background	function:Putative pheromone receptor.,miscellaneous:The gorilla and orangutan orthologous proteins do not exist, their genes are pseudogenes.,online information:No one nose - Issue 61 of August 2005,polymorphism:There seems to be a deletion in the gene coding for this protein in about 98% of the population.,similarity:Belongs to the G-protein coupled receptor 1 family.,