

VN1R4 Polyclonal Antibody

Cat No: HR1AP6863

For research use only

Overview

Product Name	VN1R4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	VN1R4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VN1R4; V1RL4; Vomeronasal type-1 receptor 4; G-protein coupled receptor GPCR27; hGPCR27; V1r-like receptor 4
Background	function:Putative pheromone receptor.,miscellaneous:The chimpanzee and gorilla orthologous proteins do not exist, their genes are pseudogenes.,online information:No one nose - Issue 61 of August 2005,polymorphism:Various VN1R4 alleles are known. The sequence shown is that of allele VN1R4*1.,similarity:Belongs to the G-protein coupled receptor 1 family.,