

## V1RL5 Polyclonal Antibody

Cat No: HR1AP6864

For research use only

### Overview

Product Name	V1RL5 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	VN1R5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VN1R5; V1RL5; Vomeronasal type-1 receptor 5; G-protein coupled receptor GPCR26; hGPCR26; V1r-like receptor 5
Background	function:Putative pheromone receptor.,miscellaneous:The chimpanzee and orangutan orthologous proteins do not exist, their genes are pseudogenes.,online information:No one nose - Issue 61 of August 2005,polymorphism:There is an allele VN1R5*2 with a truncating mutation in position 46. It was found with a frequency of about 32%.,similarity:Belongs to the G-protein coupled receptor 1 family.,