

TRH-R1 Polyclonal Antibody

Cat No: HR1AP8557

For research use only

Overview

Product Name	TRH-R1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	TRHR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRHR; Thyrotropin-releasing hormone receptor; TRH-R; Thyroliberin receptor
Background	thyrotropin releasing hormone receptor(TRHR) Homo sapiens This gene encodes a G protein-coupled receptor for thyrotropin-releasing hormone (TRH). Upon binding to TRH, this receptor activates the inositol phospholipid-calcium-protein kinase C transduction pathway. Mutations in this gene have been associated with generalized thyrotropin-releasing hormone resistance. [provided by RefSeq, Sep 2011],