

SSTR4 Polyclonal Antibody

Cat No: HR1AP8403

For research use only

Overview

Product Name	SSTR4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	SSTR4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SSTR4; Somatostatin receptor type 4; SS-4-R; SS4-R; SS4R
Background	<p>somatostatin receptor 4(SSTR4) Homo sapiens Somatostatin acts at many sites to inhibit the release of many hormones and other secretory proteins. The biologic effects of somatostatin are probably mediated by a family of G protein-coupled receptors that are expressed in a tissue-specific manner. SSTR4 is a member of the superfamily of receptors having seven transmembrane segments and is expressed in highest levels in fetal and adult brain and lung. [provided by RefSeq, Jul 2008],</p>