

## SH3TC1 Polyclonal Antibody

Cat No: HR1AP7765

For research use only

### Overview

Product Name	SH3TC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	150kDa
GeneID?Human?	SH3TC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SH3TC1; SH3 domain and tetratricopeptide repeat-containing protein 1
Background	similarity:Contains 1 SH3 domain.,similarity:Contains 9 TPR repeats.,