

## Slit3 Polyclonal Antibody

Cat No: HR1AP8347

For research use only

### Overview

|                       |  |
|-----------------------|--|
| Product Name          | Slit3 Polyclonal Antibody  |
| Source                | Rabbit   |
| Applications          | IHC-p,ELISA  |
| Species Reactivity    | Human,Mouse,Rat  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | kDa  |
| GeneID?Human?         | SLIT3  |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     | SLIT3; KIAA0814; MEGF5; SLIL2; Slit homolog 3 protein; Slit-3; Multiple epidermal growth factor-like domains protein 5; Multiple EGF-like domains protein 5  |
| Background            | slit guidance ligand 3(SLIT3) Homo sapiens The protein encoded by this gene is secreted, likely interacting with roundabout homolog receptors to effect cell migration. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012], |