

Neurexophilin-1 Polyclonal Antibody

Cat No: HR1AP6821

For research use only

Overview

Product Name	Neurexophilin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	NXPH1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NXPH1; NPH1; Nbla00697; Neurexophilin-1
Background	neurexophilin 1(NXPH1) Homo sapiens This gene is a member of the neurexophilin family and encodes a secreted protein with a variable N-terminal domain, a highly conserved, N-glycosylated central domain, a short linker region, and a cysteine-rich C-terminal domain. This protein forms a very tight complex with alpha neurexins, a group of proteins that promote adhesion between dendrites and axons. [provided by RefSeq, Jul 2008],