

# GRIN3 Polyclonal Antibody

Cat No: HR1AP6721

For research use only

## Overview

Product Name	GRIN3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GPRIN3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPRIN3; KIAA2027; G protein-regulated inducer of neurite outgrowth 3; GRIN3
Background	function:May be involved in neurite outgrowth.,