

GRIF-1 Polyclonal Antibody

Cat No: HR1AP8353

For research use only

Overview

Product Name	GRIF-1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TRAK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRAK2; ALS2CR3; KIAA0549; Trafficking kinesin-binding protein 2; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 3 protein
Background	PTM:O-glycosylated.,similarity:Contains 1 HAP1 N-terminal domain.,subunit:Interacts with GABA-A receptor and O-GlcNAc transferase (By similarity). Interacts with RHOT1/Miro-1 and RHOT2/Miro-2.,tissue specificity:Widely expressed, with highest expression in heart.,