

Formin 2 Polyclonal Antibody

Cat No: HR1AP7966

For research use only

Overview

Product Name	Formin 2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	FMN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FMN2; Formin-2
Background	<p>formin 2(FMN2) Homo sapiens This gene is a member of the formin homology protein family. The encoded protein is thought to have essential roles in organization of the actin cytoskeleton and in cell polarity. Mutations in this gene have been associated with mental retardation autosomal recessive 47 (MRT47). Alternatively spliced transcript variants have been identified. [provided by RefSeq, Mar 2015],</p>