

## Cerebellin 3 Polyclonal Antibody

Cat No: HR1AP8279

For research use only

### Overview

Product Name	Cerebellin 3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p, WB, ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID? Human?	CBLN3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CBLN3; Cerebellin-3
Background	cerebellin 3 precursor(CBLN3) Homo sapiens Members of the precerebellin family, such as CBLN3, contain a cerebellin motif (see CBLN1; MIM 600432) and a C-terminal C1q signature domain (see MIM 120550) that mediates trimeric assembly of atypical collagen complexes. However, precerebellins do not contain a collagen motif, suggesting that they are not conventional components of the extracellular matrix (Pang et al., 2000 [PubMed 10964938]). [supplied by OMIM, Aug 2009],