

## BEGAIN Polyclonal Antibody

Cat No: HR1AP8034

For research use only

### Overview

Product Name	BEGAIN Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	BEGAIN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BEGAIN; KIAA1446; Brain-enriched guanylate kinase-associated protein
Background	function:May sustain the structure of the postsynaptic density (PSD).,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subunit:Interacts with DLG4 and DLGAP1 and forms a ternary complex.,