

BM88 Polyclonal Antibody

Cat No: HR1AP7623

For research use only

Overview

Product Name	BM88 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CEND1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEND1; BM88; Cell cycle exit and neuronal differentiation protein 1; BM88 antigen
Background	cell cycle exit and neuronal differentiation 1(CEND1) Homo sapiens The protein encoded by this gene is a neuron-specific protein. The similar protein in pig enhances neuroblastoma cell differentiation in vitro and may be involved in neuronal differentiation in vivo. Multiple pseudogenes have been reported for this gene. [provided by RefSeq, Jul 2008],