

MAG Polyclonal Antibody

Cat No: HR1AP5358

For research use only

Overview

Product Name	MAG Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	MAG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MAG; GMA; Myelin-associated glycoprotein; Siglec-4a
Background	myelin associated glycoprotein(MAG) Homo sapiens The protein encoded by this gene is a type I membrane protein and member of the immunoglobulin superfamily. It is thought to be involved in the process of myelination. It is a lectin that binds to sialylated glycoconjugates and mediates certain myelin-neuron cell-cell interactions. Three alternatively spliced transcripts encoding different isoforms have been described for this gene. [provided by RefSeq, Nov 2010],