

SR-5B Polyclonal Antibody

Cat No: HR1AP5510

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | SR-5B Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IF,ELISA |
| Species Reactivity | Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 38kDa |
| GeneID?Human? | HTR5B |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of this receptor is mediated by G proteins. Probably involved in anxiety and depression. |