

Neuromedin-S Polyclonal Antibody

Cat No: HR1AP4043

For research use only

Overview

Product Name	Neuromedin-S Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	NMS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NMS; Neuromedin-S
Background	neuromedin S(NMS) Homo sapiens This gene encodes a member of the neuromedin family of neuropeptides. The encoded preproprotein is proteolytically processed to generate a biologically active neuropeptide that plays a role in the regulation of circadian rhythm, anorexigenic action, antidiuretic action, cardiovascular function and stimulation of oxytocin and vasopressin release. [provided by RefSeq, May 2016],