

Neurocan Polyclonal Antibody

Cat No: HR1AP5980

For research use only

Overview

Product Name	Neurocan Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	145kDa
GeneID?Human?	NCAN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NCAN; CSPG3; NEUR; Neurocan core protein; Chondroitin sulfate proteoglycan 3
Background	neurocan(NCAN) Homo sapiens Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration.[supplied by OMIM, Jul 2002],