

# Neurexophilin-3 Polyclonal Antibody

Cat No: HR1AP4038

For research use only

## Overview

Product Name	Neurexophilin-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	NXPH3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NXPH3; KIAA1159; NPH3; Neurexophilin-3
Background	function:May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins.,PTM:May be proteolytically processed at the boundary between the N-terminal non-conserved and the central conserved domain in neuron-like cells.,similarity:Belongs to the neurexophilin family.,tissue specificity:Highest level in brain.,