

GRIN2 Polyclonal Antibody

Cat No: HR1AP3586

For research use only

Overview

Product Name	GRIN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	GPRIN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPRIN2; KIAA0514; G protein-regulated inducer of neurite outgrowth 2; GRIN2
Background	function:May be involved in neurite outgrowth.,subunit:Interacts with activated forms of GNAO1 and GNAZ.,tissue specificity:Expressed specifically in the cerebellum.,