

Cerebellin 4 Polyclonal Antibody

Cat No: HR1AP3067

For research use only

Overview

Product Name	Cerebellin 4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	CBLN4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CBLN4; CBLNL1; Cerebellin-4; Cerebellin-like glycoprotein 1
Background	cerebellin 4 precursor(CBLN4) Homo sapiens This gene encodes a member of a family of small secreted proteins containing C1Q domains. Members of this family are involved in regulation of neurexin signalling during synapse development. The mouse homolog of the protein encoded by this gene competes with netrin to bind to the deleted in colorectal cancer receptor. [provided by RefSeq, Aug 2012],