

Cerebellin 1 Polyclonal Antibody

Cat No: HR1AP3065

For research use only

Overview

Product Name	Cerebellin 1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	CBLN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CBLN1; Cerebellin-1; Precerebellin
Background	<p>cerebellin 1 precursor(CBLN1) Homo sapiens This gene encodes a cerebellum-specific precursor protein, precerebellin, with similarity to the globular (non-collagen-like) domain of complement component C1qB. Precerebellin is processed to give rise to several derivatives, including the hexadecapeptide, cerebellin, which is highly enriched in postsynaptic structures of Purkinje cells. Cerebellin has also been found in human and rat adrenals, where it has been shown to enhance the secretory activity of this gland. [provided by RefSeq, Aug 2008],</p>