

## Cerebellin 2 Polyclonal Antibody

Cat No: HR1AP3066

For research use only

### Overview

Product Name	Cerebellin 2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	CBLN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CBLN2; Cerebellin-2
Background	similarity:Contains 1 C1q domain.,subunit:May interact with CBLN1, CBLN3 and CBLN4.,