

ALS2CR8 Polyclonal Antibody

Cat No: HR1AP2754

For research use only

Overview

Product Name	ALS2CR8 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	ALS2CR8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ALS2CR8; CARF; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 8 protein; Calcium-response factor; CaRF; Testis development protein NYD-SP24
Background	function:May be a transcription factor.,