

Neuroglycan C Polyclonal Antibody

Cat No: HR1AP4042

For research use only

Overview

Product Name	Neuroglycan C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	CSPG5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CSPG5; CALEB; NGC; Chondroitin sulfate proteoglycan 5; Acidic leucine-rich EGF-like domain-containing brain protein; Neuroglycan C
Background	chondroitin sulfate proteoglycan 5(CSPG5) Homo sapiens The protein encoded by this gene is a proteoglycan that may function as a neural growth and differentiation factor. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2011],