

DCAMKL3 Polyclonal Antibody

Cat No: HR1AP3262

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | DCAMKL3 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 73kDa |
| GeneID?Human? | DCLK3 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | DCLK3; DCAMKL3; DCDC3C; KIAA1765; Serine/threonine-protein kinase DCLK3; Doublecortin domain-containing protein 3C; Doublecortin-like and CAM kinase-like 3; Doublecortin-like kinase 3 |
| Background | catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. CaMK subfamily.,similarity:Contains 1 protein kinase domain., |