

## VGF Polyclonal Antibody

Cat No: HR1AP8608

For research use only

### Overview

Product Name	VGF Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	VGF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VGF; Neurosecretory protein VGF
Background	VGF nerve growth factor inducible(VGF) Homo sapiens This gene is specifically expressed in a subpopulation of neuroendocrine cells, and is upregulated by nerve growth factor. The structural organization of this gene is similar to that of the rat gene, and both the translated and the untranslated regions show a high degree of sequence similarity to the rat gene. The encoded secretory protein also shares similarities with the secretogranin/chromogranin family, however, its exact function is not known. [provided by RefSeq, Jul 2008].