

# Rabphilin-3A (phospho Ser237) Polyclonal Antibody

Cat No: HR1AP6380

For research use only

## Overview

Product Name	Rabphilin-3A (phospho Ser237) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	RPH3A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPH3A; KIAA0985; Rabphilin-3A; Exophilin-1
Background	rabphilin 3A(RPH3A) Homo sapiens Exocytosis of neurotransmitters and hormones is fundamental to synaptic neurotransmission and cell-cell communication. RAB3A (MIM 179390) is a small G protein that is thought to act at late stages of exocytosis, and RPH3A is a RAB3A effector (Lin et al., 2007 [PubMed 17149709]).[supplied by OMIM, Jul 2008],