

# NUMB Polyclonal Antibody

Cat No: HR1AP5161

For research use only

## Overview

Product Name	NUMB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	NUMB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NUMB; Protein numb homolog; h-Numb; Protein S171
Background	NUMB, endocytic adaptor protein(NUMB) Homo sapiens The protein encoded by this gene plays a role in the determination of cell fates during development. The encoded protein, whose degradation is induced in a proteasome-dependent manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],