

M6A Polyclonal Antibody

Cat No: HR1AP4885

For research use only

Overview

Product Name	M6A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	GPM6A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPM6A; M6A; Neuronal membrane glycoprotein M6-a; M6a
Background	similarity:Belongs to the myelin proteolipid protein family.,