

# Aladin Polyclonal Antibody

Cat No: HR1AP4883

For research use only

## Overview

Product Name	Aladin Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	AAAS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AAAS; ADRACALA; GL003; Aladin; Adracalin
Background	<p>aladin WD repeat nucleoporin(AAAS) Homo sapiens The protein encoded by this gene is a member of the WD-repeat family of regulatory proteins and may be involved in normal development of the peripheral and central nervous system. The encoded protein is part of the nuclear pore complex and is anchored there by NDC1. Defects in this gene are a cause of achalasia-addisonianism-alacrima syndrome (AAAS), also called triple-A syndrome or Allgrove syndrome. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010],</p>